

SUGGESTIONS

Albert y Ferran Adrià

Texturas **Suggestions**

Reverse blackcurrant Spherical



For the Algin bath:

- 1500 g mineral water
- 7.5 g Algin

- 1 Dissolve the Algin in the water with a Túrmix and leave to rest in the refrigerator.

For the spherical mixture:

- 100 g water
- 9 g Gluco
- 0.9 g Xantana
- 200 g Søbogaard blackcurrant nectar

- 1 Mix the water and the Gluco, add the Xantana with a Túrmix.
- 2 Leave the mixture to lose the air and add the blackcurrant nectar.
- 3 Make the sphericals in the Algin bath and leave to cook for a minute and a half.
- 4 Take out of the bath and wash in water.
- 5 Once washed, keep them in blackcurrant nectar.

Others

- c.s. Søbogaard blackcurrant nectar
- 1 set of Eines

- * this spherical goes perfectly with cocktails as in this case a Triple Sec liqueur.